

ParticleMaster

Intelligent Imaging
for Particle & Droplet Sizing




LAVISION

FOCUS ON IMAGING



ParticleMaster Shadow

ParticleMaster Shadow is based on backlight illumination and high-magnification imaging. The shadow of particles in the focal plane of the optics is imaged.




detector: long distance microscope with high resolution CCD


particles, droplets

illumination optics

light source: (pulsed) laser or (flash) lamp



air bubbles in water



droplets and ligaments in air

Information

- ▶ particle size (d)
- ▶ particle position (x, y)
- ▶ particle shape (excentricity)
- ▶ statistics, histograms (D_{10} , D_{32} , D_{v50})
- ▶ velocity (v_x , v_y)
- ▶ density
- ▶ mass flux

Applications

- ▶ liquid sprays (water, fuel, paint, emulsions)
- ▶ spray breakup (ligaments, breakup region)
- ▶ powder, solid particles (alloys, ceramics)
- ▶ bubbles (heat exchangers, industrial processes)


System Components

- ▶ ParticleMaster Shadow Sizing software
- ▶ CCD camera*
- ▶ long distance microscope or macro lens
- ▶ laser*, laser diode flasher or flash lights, illumination optics

*standard PIV components

Upgrades

- ▶ PIV or planar LIF systems for multi-phase flows



particle size histogram



ParticleMaster IMI

ParticleMaster IMI is based on **Interferometric Mie Imaging**. The spatial Mie scattering intensity distribution is recorded. Size information is obtained from defocused imaging.



$$d = \frac{2 \cdot \lambda}{\Delta\phi} \cdot \left[\cos\left(\frac{\phi}{2}\right) + \frac{n \cdot \sin\left(\frac{\phi}{2}\right)}{\sqrt{1 + n^2 - 2 \cdot n \cdot \cos\left(\frac{\phi}{2}\right)}} \right]^{-1}$$


Reference: W. J. Glantschnig, S.-H. Chen, Appl. Opt., 20 (14), 1981.

Information

- ▶ **particle size** (d)
- ▶ **particle position** (x, y, z)
- ▶ **velocity** (v_x, v_y, v_z)
- ▶ **statistics, histograms** (D₁₀, D₃₂, D_{v50})
- ▶ **density**
- ▶ **mass flux**

Applications

- ▶ transparent sprays (water, fuel, pharma-sprays)
- ▶ droplet clouds (evaporation and condensation)
- ▶ bubbles (heat exchangers, industrial processes)



Reference: A. Graßmann, F. Peters, personal communication, 10th Fachtagung GALA, Rostock, 2002, 11th Fachtagung GALA, Braunschweig, 2003.

Special IMI Features

- ▶ auto-detection and fringe pattern analysis (DaVis IMI Sizing)
- ▶ FFT with sub-pixel accuracy & advanced evaluation algorithms
- ▶ centroiding algorithm for high-accuracy positioning and velocity determination
- ▶ high-aperture macro lens


System Components

- ▶ **ParticleMaster IMI Sizing** software
- ▶ CCD camera*, macro lens
- ▶ laser*, sheet optics*


*standard PIV components

Upgrades

- ▶ PIV or planar LIF systems for multi-phase flows



ParticleMaster Application Matrix



	Shadow	IMI
Particle size in μm	> 5	2 - 200
Dynamic range	1:100	1:20
Typ. field of view in mm^2	> 1	50 - 500
Typ. working distance in mm	50 - 1500	50 - 300
Spray density n	$10 n_{\text{IMI}}$	n_{IMI}
Opaque particles	yes	no
Shape	yes	no
Velocity	yes	yes
Morphology	yes	no

LaVisionUK Ltd

Downsview House / Grove Technology Park
Grove / Oxon / OX12 9FF / United Kingdom
E-Mail: sales@lavidionuk.com
www.lavidionUK.com
Phone: +44-(0)-870-997-6532
Fax: +44-(0)-870-762-6252

LaVision GmbH

Anna-Vandenhoeck-Ring 19
D-37081 Goettingen / Germany
E-Mail: info@lavidion.com
www.lavidion.com
Tel.: +49-(0)5 51-9004-0
Fax: +49-(0)551-9004-100

LaVision Inc.

211W. Michigan Ave. / Suite 100
Ypsilanti, MI 48197 / USA
E-Mail: sales@lavidioninc.com
www.lavidion.com
Phone: (734) 485 - 0913
Fax: (240) 465 - 4306