

Bioaerosol Sampling

High efficiency detection of
airborne virus using a
condensation growth
sampling method

Evidence for airborne transmission of SARS-CoV-2



Pat Keady, President



Aerosol Devices Inc

pkeady@aerosoldevices.com

The Problem



Highly Mobile Society



“There is no
evidence of
aerosol
transmission”

WHO still will not recognize
aerosols as a significant route of
the SARS-CoV-2 virus
transmission

When scientists say “there is no evidence”,
sometimes it’s because that evidence is
genuinely hard to collect.

Anecdotal evidence: rapid disease spread



Diamond-Princess cruise ship




Church choir practice



San Quentin Prison



Minnesota bars



Guangzhou China restaurant





Anecdotal evidence: virus found on air vents, in nurse changing rooms, and in hallways outside of infected patient rooms.

Added to the
WHO website
just yesterday,
July 9th !

There have been reported outbreaks of COVID-19 reported in some closed settings, such as restaurants, nightclubs, places of worship or places of work where people may be shouting, talking, or singing. In these outbreaks, **aerosol transmission**, particularly in these indoor locations where there are crowded and inadequately ventilated spaces where infected persons spend long periods of time with others, **cannot be ruled out. More studies are urgently needed** to investigate such instances and assess their significance for transmission of COVID-19.

<https://www.who.int/news-room/q-a-detail/q-a-how-is-covid-19-transmitted>

What is an aerosol?



Vincent, J.H. J Env Monitoring, <https://doi.org/10.1039/B509617K>



Do we have
a common
language?


Seeing a bright side to the COVID-19 pandemic,
one aerosol scientist said to another:

“At least now people won’t think ***aerosols*** are
just spray cans and graffiti.”

“Yeah, now they think aerosols are ***viruses!***”

Pathogen emissions expelled from human respiratory tract

- Coughing and Sneezing
 - Singing
 - Loud talking
 - Laughing
 - Just breathing
-
- Particles expelled are a continuum of sizes





Super spreaders

Some people, known as “super spreaders”, happen to be especially good at exhaling fine material, producing 1,000 times more than others.

Pre-symptomatic and asymptomatic cases can be considered “stealth” spreaders.

Bacteria from respiratory droplets spread by...

...one sneeze



...singing
(1 min)



... talking
(1 min)



... two
coughs



No Mask

Masked

Dr. Richard Davis, Providence Sacred Heart Medical Center



Note: In this demonstration presence of bacteria (not viruses) on plates is only meant to be a proxy for microbes present in respiratory droplets.

Likely, smaller aerosolized droplets (that could carry viruses like SARS-CoV-2) are also produced by coughing, sneezing etc. and that these would travel further and stay in the air longer than larger respiratory droplets.

Masks Work

Which masks work best?




- **Best:** Powered air purifying respirator (PAPR) - what fire fighters use
- **Better:** Elastomeric face masks with cartridge filters – high efficiency and great face seal, fit tested - for dangerous occupational exposures
- **Good:** N95 (FFP2) Filtering face pieces, fit tested – what medical doctors/nurses use
- **Fair:** Surgical masks and home-made cloth masks – remove large expelled droplets and are most protective of others, but can still offer some protection to the wearer if worn correctly
- **Poor:** NO MASK

Arguably the best description for how a N95 (or FFP2) mask works – accurate AND understandable: <https://youtu.be/eAdanPfQdCA>

COVID-19 Airborne Transmission Tool Available

New model estimates
COVID-19
transmission in
classrooms, buses,
protests, more



Developed by J Jimenez, Univ Colorado
Boulder, with input from many aerosol
scientists

<https://cires.colorado.edu/news/covid-19-airborne-transmission-tool-available>

<https://docs.google.com/spreadsheets/d/16K1OQkLD4BjgBdO8ePj6ytf-RpPMIJ6aXFg3PrIQBbQ/edit#gid=519189277>

Missing evidence: what we still don't know about airborne SARS-CoV-2

- Particle size distribution and viral content from exhaled breath
- Viability of airborne viruses
- Airborne transportation and exposure
- Dose response – how many viable virions cause infection?
- Viral shedding load at various stages of infection

Knowledge will lead to more effective non-pharmaceutical interventions

Engineering Controls

- Ventilation
- Physical barriers
- Filtration
- UV light
- Germicidal sprays

Social distancing

- Is 2m distance far enough? Is 1m sufficient?

Face Masks

- N95 vs surgical masks
- Homemade cloth coverings
- Fit Testing




We are blind to the
biohazards
in the air we breathe

Bioaerosols are Difficult to Sample

- Size-dependent collection efficiency
- Reduced collection into liquids for hydrophobic organisms
- Dilute liquid samples
- Prone to contamination
- No regulatory standard method

Common Sampling Methods



Wet-walled
cyclone



Filter
Sampling
System



Sampling Viable Bioaerosols

Viable organisms inactivated by the sampling method

- Dessication
- Mechanical stress
- Thermal shock



A New Approach for Sampling Aerosols for Biological Analysis



SPOT SAMPLER™
Aerosol Particle Collector

- ✓ High collection efficiency
- ✓ Concentrated sample
- ✓ Time-resolved sampling
- ✓ Maintains viability
- ✓ Instant genomic-DNA/RNA preservation

BioSpot-VIVAS™
Bioaerosol Sampler



Patented technology with exclusive license from Aerosol Dynamics Inc.

Brief History of Water-based Condensation Particle Growth

First systems: mix steam into airstream, then cool


Maze Collector (Simon and Dasgupta, 1995)

Steam-jet Collector (Khlystov et al, 1995)


Particle-in-Liquid Sampler (Weber et al, 2001)

These methods subject the sample to high temperatures

2003, S. Hering (Aerosol Dynamics Inc.) introduced technology that allows particle growth through laminar flow water condensation





New Design



Condensation Growth Tube = CGT

A. Eiguren-Fernandez et al., presented at RICTA 2015, the 3rd Iberian Meeting on Aerosol Science and Technology in Elche, Spain, June 29-July 1, 2015.
Susanne V. Hering , Steven R. Spielman & Gregory S. Lewis, Aerosol Science and Technology, 48:4, 401-408: 2014.

Condensation Growth Tube (CGT) Capture




Moderate sample flow
temperatures never exceed 30°C.
Exit flow temperature <18°C;
dewpoint < 20°C.

Supersaturation levels of 120-
140% activate condensation
growth on particles as small as
5 nm.


Droplets grown to nominal 3μm
diameter are easily captured by
bounce-free, soft inertial
impaction on to a solid surface
or into liquid

Growth Tube - Efficient over all Particle Sizes



Data for wet-walled cyclone from McFarland et al (2010); BioSampler data from Hogan et al (2005) and Willeke et al. (1998).


Collection of MS2 virus aerosol




Condensation Growth Tube = CGT

M. Pan et al., Journal of Applied Microbiology; 120, 805-815: 2016.


Collection of infectious *viable* influenza H1N1 virus aerosol




Sampling period (min)	Number of tests	Collection efficiency of infectious H1N1	
		BioSampler	CGT
5	7	$4.4 \pm 2.6\%$	$67 \pm 10\%$
10	3	$6.0 \pm 0.8\%$	$83 \pm 12\%$
15	3	$8.2 \pm 4.3\%$	$80 \pm 07\%$
All Runs	13	$5.6 \pm 3.0\%$	$74 \pm 12\%$



Collection of airborne viable bacteria and yeast





	BioSampler (CFU/L)	CGT (CFU/L)
YTG	147 ± 20	471 ± 175
PBS	146 ± 7	415 ± 61






Condensation Capture into Genomic Preservative



Condensation Growth Tube = CGT
No Condensation Growth Tube = NCGT

M. Nieto Caballero et al., 2019, High Fidelity Recovery of Airborne Microbial Genetic Materials by Direct Condensation Capture into Genomic Preservatives.

Bioaerosol Sampling and Analysis



Biological Assays

Fluorescence spectroscopy

Microscopy techniques

qPCR, DNA/RNA Sequencing

Culture

Microfluidic Sensing

MALDI-TOF spectrometry

Raman spectroscopy



FINNISH
METEOROLOGICAL
INSTITUTE



RTI
INTERNATIONAL

UF | UNIVERSITY of
FLORIDA



Washington
University
in St. Louis



Mount
Sinai



UNIVERSITÉ
LAVAL



Early Adopters





U.S Food and Drug Administration



Future Developments



IH/IAQ Consult
Health care
Public health
Agribusiness
Eldercare Homes
Transportation
Restaurants/Bars



Low cost, simple sampler for genomic analysis – qPCR, RNA/DNA sequencing



Sampler integrated with an electrochemical biosensor targeted for a specific virus or bacteria; focus is initially on SARS-CoV-2

<https://www.youtube.com/watch?v=ujFhI42R2-8>

Condensation Growth Tube Capture – Advantages for Bioaerosols

- Gentle capture into liquid maintains viability
- Capture into a genomic stabilizer instantly preserves DNA/RNA
- Choice of collection terminus – liquid or solid substrate
- High, uniform collection efficiency independent of particle size, shape, composition, or hydrophobicity
- Wide particle size range from 5nm to 10 μ m
- Concentrated sample
- No temperature extremes – mimics the environmental conditions in the human lung
- No particle bounce or re-aerosolization
- Minimal sample handling → low contamination potential

For a
healthier
future,
together.

Thank you!



Aerosol Devices Inc

Pat Keady

pkeady@aerosoldevices.com

970-631-5950