


BRECHTEL

**Solutions for your
research challenges**

ISO Isokinetic Inlet System

Model 1200




When you need to get a representative aerosol sample to your instruments on an aircraft or mobile platform

Features:

- Fully automated, hands-off operation
- Transmission efficiency >90% for particle aerodynamic diameters <10 μm
- Wind tunnel tested and validated
- Inlet sampling manifold includes instrument sample pickoffs
- Custom mounting solutions that account for aircraft angle of attack (optional)
- Aerodynamic, twin-diffuser prevents turbulent boundary layer separation to minimize particle losses
- Low power blower to draw and actively control 150 lpm sample flow (up to 100 lpm of sample air flow for instruments)
- Anti-icing system

**www.brechtel.com
sales@brechtel.com**

Providing Aerosol Measurement Solutions



Specifications

Parameter	Value
Particle aerodynamic diameter size range	0.005-10 μm
Sample flow at diffuser tip (@100 m/s)	300 lpm
Sample flow to cabin	150 lpm
Maximum sample flow to instruments	100 lpm
Control and DAQ frequency	1 Hz
Anti-icing power (supplied by platform)	900 watts @ 28 VDC
Chassis power (supplied by platform)	56 watts @ 28 VDC
Rack-mountable electronics chassis size	19 x 6.5 x 12 in 48.25 x 16.5 x 30.5 cm
Electronics chassis weight	20 lb/9 kg
Total system weight	70 lb/31.75kg
Operating temperature range	-40 to 45 °C
Operating pressure range (absolute)	200-1,000 mb


Copyright ©2020 specifications are subject to change without notice. BMI assumes no responsibility for inaccuracies in this document or for any obligation to update information in this document. BMI reserves the right to change, modify, transfer or otherwise revise this publication without prior notice.

*Some products may be shown with optional accessories, which are sold separately. Items shown may not be to scale.

Publication:

R. Zaveri, C. M. Berkowitz, F. J. Brechtel, M. K. Gilles, J. M. Hubbe, J. T. Jayne, L. I. Kleinman, A. Laskin, S. Madronich, T. B. Onasch, M. S. Pekour, S. R. Springston, J. A. Thornton, A. V. Tivanski, and D. R. Worsnop (2009). [Nighttime chemical evolution of aerosol and trace gases in a power plant plume](#), J. Geophys. Res., 115, D12304, doi:10.1029/2009JD013250.

Schematic of the ISO Inlet



ISOKINETIC INLET ASSEMBLY SCHEMATIC

Applications

- Plume studies of power plants
- Megacity plume studies
- Volcanic ash studies
- Below cloud aerosol studies for cloud condensation nucleus research
- Visibility reduction studies
- Climate change research
- Vertical profiling of aerosols
- Air quality research
- Mobile laboratories

How to Order

Part No.	Description
1200	Isokinetic Sampling (ISO) Inlet System
ISO-SP025	Sample pick-off port with 1/4" OD tube
ISO-SP0375	Sample pick-off port with 3/8" OD tube
ISO-SP05	Sample pick-off port with 1/2" OD tube
ISO-90	Ninety degree bend for 1.5" OD main sample manifold line at mounting plate
ISO-Kit	Maintenance Kit for 1200 ISO